
MACRO-LOCK® POST X-RO® ILLUSION®

- Ⓒ GB INSTRUCTION FOR USE
- Ⓒ FR MODE D'EMPLOI
- Ⓒ BR INSTRUÇÕES DE USO
- Ⓒ PT INSTRUÇÕES DE USO
- Ⓒ DE GEBRAUCHSANWEISUNG
- Ⓒ ES INSTRUCCIONES DE USO
- Ⓒ GR ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
- Ⓒ 中文 使用说明
- Ⓒ KR 사용설명서
- Ⓒ RU ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ

NOT4590XX00-201411-04

INSTRUCTIONS FOR USE

DESCRIPTION

MACRO-LOCK POST® ILLUSION® X-RO® is a serrated taper post, length is 17,5 mm.

MACRO-LOCK POST® ILLUSION® X-RO® is a yellow translucent fiber post with a colored resin matrix which becomes translucent in the mouth. The post color will re-appear on command under cold air, water spray if re-treatment is required.

The radiopacity level of the **MACRO-LOCK POST® ILLUSION® X-RO®** is 340 % Al equ. (ISO 4049)

	● 1	● 2	● 3	● 4	● 5	● 6
Ø apical tip	0.80	0.80	1.00	1.00	1.30	1.30
Ø Post head	1.35	1.47	1.67	1.83	2.05	2.20
Starter drill	#1	#1	#2	#2	#2	#2
Finishing Drill	#1	#2	#3	#4	#5	#6

INDICATIONS FOR USE

In case of insufficient residual tooth structure (<4 mm), the post is needed to support and secure the coronal restoration.

CONTRAINDICATIONS

If coronal tooth structure is less than 1.5 mm.

PRECAUTIONS

The post must be cleaned with alcohol before insertion in the canal. Avoid touching the posts with your fingers after cleaning.

Shortening the post should be done outside of the mouth.

The use of rubber dam is recommended.

The crown must cover the preparation of at least 1.5 mm of healthy dentin in order to get the ferrule effect.

ADVERSE REACTIONS

None known to date.

STEP-BY-STEP INSTRUCTIONS

1. Root canal preparation: remove gutta-percha to preplanned depth with a Gates-Glidden, Peeso reamer or Largo. Radiographic verification is recommended.
2. Select the correct MACRO-LOCK POST® **ILLUSION® X-RO®** size. Determine the post size according to the anatomical situation of the tooth using the radiograph and the plastic calibrating card.
3. Determine the length of the coronal build-up, which will account for 1/3 of the final post length. Remove the root filling with the Macro-Lock Starter Drill selected (rotation speed 800-1200 rpm) as deep as necessary for the post to be inserted to 2/3 of the root length. At least 4 mm of root canal filling should remain in the apical region. In case of curved canal, this calculated length is not suitable.
4. Shape the canal with the Finishing Drill corresponding to the selected **MACRO-LOCK POST® ILLUSION® X-RO®** size, (rotation speed 1000-2000 rpm).
5. Try in the post, to check that the post is properly seated in the canal
6. Shorten the post to its final length with a diamond disc outside of the mouth. **Never use a crimping instrument, such as wire cutters, as the pressure can destroy the structure of the post.**
7. Clean the post with alcohol.
8. **Apply etchant in the post space and to the exposed dentin for 15 seconds*. Rinse for 10 seconds. Remove the excess water with paper points but leave surface moist. Do not put the etchant in contact with the gingiva.**
9. Apply 2 coats* of the primer with a thin brush (e.g. Composibrush RTD) in the post space. Remove excess using paper points and gently air dry all surfaces. Enamel and dentin surfaces should have a uniform, glossy appearance. If not, repeat application. Light-cure the primer* inside the post space for 10 to 20 seconds.
10. Apply a single coat of primer to the post. Gently air dry 5 seconds and light-cure the post for 10 to 20 seconds (if needed) outside the mouth*.
- 11a. If using a dual cure resin cement for both cementation and core build up such as Corecem™ (recommended technique): Mix the cement and apply it in the post space from bottom to the top using the root canal tips. Seat the post immediately. Light-cure for 40 to 60 seconds. Continue applying the core build up using a plastic Coreform™ matrix, or directly on the tooth.
- 11b. If using separate cement and a core build up resins : Mix the cement* and apply it to the post and then in the post space using a Lentulo spiral or a syringe tip. Seat the post immediately. Remove excess cement with appropriate instruments. If using dual-cure cement, light-cure for 40 to 60 seconds, applying gentle pressure on the post with the tip of the curing light probe. Apply 2 coats of bonding agent* to the exposed post(s), cement and involved tooth structure. Remove excess and air-dry with jet of air. Light-cure for 20 seconds (if needed).

- 11c.i. If using a self etching cement*, there is no need to etch and prime the post space. Mix the cement and apply it in the post space from bottom to the top using a root canal tip. Seat the post immediately. Remove excess cement. Light cure if needed.
- 11c.ii. After etching, rinsing, drying of involved tooth structure, apply 2 coats of bonding agent* to the exposed post(s) and coronal part. Remove excess and air dry with jet of air. Light cure for 20 seconds (if needed).
12. Directly model the build-up using a composite core build-up material*. R.T.D. recommends that the coronal end of the post be covered by core composite.

OTHER RECOMMENDATIONS

A thick rubber dam can momentarily prevent the post to change color.

MACRO-LOCK POST® ILLUSION® X-RO® are single use devices. They must not be re-used to avoid the risk of contamination.

MACRO-LOCK POST® ILLUSION® X-RO® can be cleaned with alcohol. They can neither be sterilised nor disinfected in thermoisinfecteur.

The use of autoclave is not recommended.

A re-access carbide kit is available if needed. Contact your RTD distributor.

Starter Drill and Finishing Drill: these instruments are supplied non-sterile. Disinfect and sterilize instruments before each use. Disinfect instruments with disinfecting or cleaning agents adapted for rotary instruments*. Do not disinfect instruments in thermoisinfecteur. Inspect instruments visually before each use: change to new drills when they are damaged, corroded or dull. The universal drill may be used only to remove root canal filling partially, and the finishing drills may be used only to prepare root canal.

Drills must be sterilized individually in an autoclave with the following settings**:

- Wrap individually in disposable sterilization pouches or tubing that comply with the regulation standard in force in the country of use (e.g EN ISO 11607-1),
- Autoclave: type B complying with regulation standard in force (e.g EN 13060),
- **Sterilization temperature: 134°C - sterilization time: 18 minutes.**

Use a new drill after 12 -15 uses.

Store the sterilized components in a dry, dust-free place.

If the packaging's integrity appears to be compromised, before using again, place the component in a new pouch and re-sterilize according to the protocol defined before.

For profoundly flared or ovoid canals, the use of accessory post, such as **Fibercone®**, is indicated for best results.

SAFETY AND LIABILITY

WARRANTY

RTD warrants this product will be free from defects in material and manufacture. RTD makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and RTD's sole obligation shall be repair or replacement of the RTD product.

LIMITATION OF LIABILITY

Except where prohibited by law, RTD will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only.

** Refer to the manufacturer's instructions for use.*

*** Please refer to your autoclave manufacturer for specific sterilization instructions if needed*

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