CORECEM®

CORECEM® ILLUSION®

- Ⓜ GB INSTRUCTION FOR USE
- Ⓜ FR MODE D'EMPLOI
- Ⓜ DE GEBRAUCHSANWEISUNG
- Ⓜ GR ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
- Ⓜ ES INSTRUCCIONES DE USO
- Ⓜ RU ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
- Ⓜ PT INSTRUÇÕES DE USO
- Ⓜ BR INSTRUÇÕES DE USO

INSTRUCTIONS FOR USE

DESCRIPTION

CORECEM® has been designed to be used both for cementation of RTD fiber posts and building of the Core Build up to obtain quickly and easily a monobloc durable restoration.

CORECEM® is a dual-cured flowable nano hybrid composite.

Delivered in auto-mixing 1:1 MINIMIX (5ml) or AUTOMIX (25ml)-syringes it can be easily dispensed and applied directly into the root canal and in the COREFORM™ matrix by using root canal tips. **CORECEM®** exhibits a short setting time without high heat generation. The dual cure properties enables the dentist, also to carry out cementations and core build ups in cases where a light cure cannot be guaranteed to be sufficient.

SHADE

CORECEM® is available in two shades: A1** and **ILLUSION Blue**. Except the shades, they have similar characteristics

CORECEM® meets the requirements of the international standard ISO 4049.

INDICATIONS

- Fiber Posts cementation,
- Core build up.
- **ILLUSION Blue** is for posterior restoration, its' contrasting blue color disappears at body temperature.

CONTRAINDICATIONS

- Insufficient residual dentin : at least 1,5 mm of tooth structure round the preparation.

PREPARING THE MIXING SYRINGES

- Remove the cap of the **CORECEM®**-syringe and throw it away (do not use it again!). It is replaced by a special 1:1 mixing cannula. Turn the cannula 90° until it locks in position. The material is now ready for application.
- Note: Discard the first 2-3 mm (about the size of a pea) of the extruded material. This must be done for each new mix.

- Store used **CORECEM®**-syringe with fixed used mixing cannula in the dark.
- Do not pull back on the plunger after dispensing. This may cause plugging of the syringe.

Working time in the self-cured mode (23°C (74°F)): 1:30 minutes from start of mixing.

Setting time in the self-cured mode (23°C (74°F)): 5:30 minutes from start of mixing.

POST CEMENTATION

- Prepare the post space according to manufacturer's instructions. Etching and Primer are required.
- Dispense a small amount 2-3mm to eliminate the previous mixing and voids.
- Apply **CORECEM®** inside the root canal using the root canal tip (recommended). Filling should occur from bottom upward to prevent voids.
- Seat the post, light-cure for at least 20s. (400mW/cm²)

CORE BUILD-UP

- Fill in **CORECEM®** into the COREFORM™ matrix or fill in directly into the preparation. Filling should occur from bottom upward to prevent voids
- Light-cure for at least 40s 400mw/cm² or allow the system to self-cure for at least 3 minutes
- Remove the matrix when the material is completely set.
- Start the final crown preparation.
- For **CORECEM® ILLUSION** use cold water spray during trimming. After finishing, blue color will become dentin shade like A1.

OTHER RECOMMENDATIONS

- Do not use any resin to adjust viscosity of composite restorative material.
- Do not use in combination with substances containing eugenol because eugenol can impair the hardening of the composite and causes discoloration. Do neither store the composite material in proximity of eugenol containing products, nor let the composite allow to come into contact with materials containing eugenol.
- Contact of resin pastes with skin and gingival tissue should be avoided, especially by anyone having known resin allergies to methacrylates.
- Use a technique isolating as a rubber dam or matrix band...
- Store unopened material at a dark and cold place. Temperature should not exceed 20 °C (68 °F). Opened cartridges should be used up within 3 months. The shelf-life is two years. Do not use after expiration date.

SAFETY AND LIABILITY

WARRANTY

RTD warrants this product will be free from defects in material and manufacture. RTD makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and RTD's sole obligation shall be repair or replacement of the RTD product.

LIMITATION OF LIABILITY

Except where prohibited by law, RTD will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only.

MANUFACTURER

RTD

3 rue Louis Neel
38120 St Egrève
France

www.rtd.fr
rtd@rtd.fr

**** VITA is a registered trademark of VITA ZAHNFABRIK**

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