

# MATCHPOST™



Translucent glass fiber

## ALL THESE POPULAR FEATURES


The best value for your money!

- 
- TAPERED TO CONSERVE DENTINE
  - RADIOPAQUE
  - FOUR SIZES. TO FIT EVERY TOOTH
  - TRANSLUCENT
  - STRONG
  - DURABLE

- Radiopaque (2,00mm Al (ISO 4049)
- Elastic Modulus similar to that of dentin, to prevent root fractures.
- Micro-retentive surface, for high bond strength.
- High strength, for clinical durability.
- Biocompatible and cannot corrode.
- Aesthetic.




MATCHPOST Radiopacity  
(Courtesy Dr. B. Duret)



**COMPOSITION:**

Glass fiber: >60%  
Epoxy resin matrix: <40%

**Clinical case** (Courtesy Dr Bertoldi)




**4800001** 12 Posts Ø 1.0 - 1.2  
+ Drills




**4820210** 1 Drill Ø 1.0



**4810010**  
10 MATCHPOST Ø 1.0



**SUGGESTED POST SIZE**



3 rue Louis Neel Technoparc Espace Gavanière  
38 120 ST EGREVE - FRANCE  
Tél.: + 33 (0)4 76 56 56 66 Fax: + 33 (0)4 76 56 11 77  
E-mail: rtd@rtd.fr / www.rtd.fr



0459  
2006